

2011 Military Health System Conference

Healthcare Quality and Patient Safety Innovations: Lessons from the Field

Improving the High-Level Disinfection Process of Vaginal Ultrasound
Probes

The Quadruple Aim: Working Together, Achieving Success

LCDR Wendy Anne Cook, NC, USN

25 January 2011



Naval Hospital Bremerton,
Washington

Project Background



- Proactive evaluation of high-level disinfection (HLD) process
- Vaginal ultrasound probes selected
- 5 clinical areas
- Lean Six Sigma and H-FMEA methodologies
- Multidisciplinary team
- 6 months



Key Findings



- HLD is a complicated process
- Multiple guidelines to follow
- Measurement tool needed
- Obsolete equipment in use
- Centralization of HLD not practical
- Staff apprehensive to fully immerse probes
- Practical, economical solutions identified
- Leadership support imperative




High-Level Disinfection Audit Tool



Note: Design and Process sections omitted to improve visibility

STANDARD	POINT S	YES	NO
Equipment			
5. Is a clock or timer in the room?	1		
6. Are gloves in the room?	1		
7. Are goggles or eye protection in the room?	1		
8. Are gowns in the room?	1		
9. Is the disinfectant solution container the appropriate size?	1		
10. Is there a rinse container?	1		
11. Is the rinse container the appropriate size for high volume rinse (2 gallons)?	1		
12. Are step-by-step instructions clearly posted on the wall?	1		
13. Are equipment manufacturer guidelines readily available?	1		
14. Are disinfectant solution manufacturer guidelines readily available?	1		
Overall			
25. Are the equipment manufacturer guidelines followed?	2		

Compliance Based on High-Level Disinfection Audit Tool



Failure Modes and Effects Analysis




Failure Modes and Effects (Risk)		Baseline Risk				Risk Reduction	Revised Risk			
Process Step	Failure Mode & Effect	Sev	Occ	Det	RPN	Recommended Actions	Sev	Occ	Det	RPN
Disinfectant tested once daily in AM	Disinfectant not tested prior to each use, concentration ineffective	10	5	10	500	Test disinfectant prior to each use	10	1	10	100
3/4 of probe handle immersed in disinfectant	Probe partially immersed in disinfectant, contaminated area of probe not disinfected	5	10	5	250	Fully immerse probe handle and 12-18 inches of cord in disinfectant	5	1	5	25
Probe air-dried following disinfection	Bacterial growth encouraged due to delayed drying time, probe contaminated	5	5	10	250	Dry probe with sterile gauze following disinfection	5	1	10	50
Probe returned to exam room uncovered	Uncovered probe transported through clinic, probe contaminated	1	5	10	50	Cover probe with clean cloth for transport and storage	1	1	10	10
Total					1050					185

Results



Compliance Based on High-Level Disinfection Audit Tool



Post Improvement



Obsolete wall unit removed,
wall repaired

Test log

Posted
instructions

High-level disinfectant
containers (allow full
immersion of probe and
12-18 inches of cord,
probe soaks for 12-60
minutes)

High volume rinse
containers
(~2 gallons,
probe soaks in
fresh water rinse
for 1 minute x3)

Table added
for counter space

Utility sink



Benefits



- Improved compliance
- Decreased process variation
- Safer patient care environment
- Easily replicated in other clinical areas
- Successful use of LSS to improve patient care
- Team satisfaction
- Easily sustained

